












Main Menu File Experiment Types Setup Reports About

PRESARIO-5834 02-04-2001 04:07:00 AM

Main Menu Options

 Experiment Setup	 Data Upload
 Experiment Sessions	 Data Backup
 Calibration	 Security
 Synchronize Databases	 Event Log
 Lookup Settings	 Alarm Log
	 Search

Left Column Callouts:

- Select Experiment (points to Experiment Setup)
- View sessions management (points to Experiment Sessions)
- Calibrate instrumentation (points to Calibration)
- Configure Data Synchronization (points to Synchronize Databases)
- Configure Experiment (points to Lookup Settings)

Right Column Callouts:

- Upload Data to corporate wide Discovery DB (points to Data Upload)
- Backup Setup and configuration (points to Data Backup)
- Admin level Security access and configuration (points to Security)
- View Event Log (points to Event Log)
- View Alarm Log (points to Alarm Log)
- Advanced Search and Data mining reporting. (points to Search)

Navigation buttons: [Back] [Home] [Forward] [Refresh]